

PET2BIDS: a library for converting Positron Emission Tomography data to BIDS


Anthony Galassi ¹, **Martin Norgaard** ^{2,3,4}, **Adam G. Thomas** ¹, **Gabriel Gonzalez-Escamilla** ⁵, **Claus Svarer** ², **Chris Rorden** ⁶, **Granville J. Matheson** ^{7,8,9}, **Gitte M. Knudsen** ², **Robert B. Innis** ¹, **Melanie Ganz** ^{2,3}, **Cyrus Eierud** ¹⁰, **Murat Bilgel** ¹¹, and **Cyril Pernet** ²

¹ National Institutes of Health, Bethesda, MD, United States ² Neurobiology Research Unit, Rigshospitalet, Copenhagen, Denmark ³ Department of Computer Science, University of Copenhagen, Copenhagen, Denmark ⁴ Department of Psychology, Stanford University, CA, United States ⁵ University Medical Center of the Johannes Gutenberg University Mainz, Mainz, Germany ⁶ Department of Psychology, University of South Carolina, Columbia, SC, United States ⁷ Mailman school of Public Health, Columbia University, New York, NY, United States ⁸ Department of Clinical Neuroscience, Karolinska Institutet and Stockholm County Council, Stockholm, Sweden ⁹ Karolinska Institutet ¹⁰ TRenDS Center, Georgia State University, Atlanta, GA, United States ¹¹ National Institute on Aging Intramural Research Program, Baltimore, MD, United States

Features

- Matlab and Python
- DICOM and ecats 7+
- Metadata handling
- Blood data conversion

→ the library converts files (.nii .json) and names files following BIDS, you still have to structure your data (folders and the like)



Update existing json

- *updatejsonpetfile.m or update_json_pet_file.py*
- Internal routines that check and possibly add tracer info (e.g. are all the values consistent with expected units Bq, MBq, GBq, g, ug, mol, umol)
- Validates required, recommended and optional fields
- Apply proper names to 'identified' reconstruction method (heuristics)

Update existing json from excel

- Convert preformatted xls, xlsx, csv, tsv

participant_id	InjectedRadioactivity	SpecificRadioactivity	InjectedMass	TimeZero	DecayCorrectionFactor
D:\BIDS\ONP\BIDS-converter\spreadsheet_conversion\many_subjects_sheet\sub-01\pet	397	59,31	1,896978587	13.42.15	
D:\BIDS\ONP\BIDS-converter\spreadsheet_conversion\many_subjects_sheet\sub-02\pet	551	25,28	6,176954114	09.12.12	
D:\BIDS\ONP\BIDS-converter\spreadsheet_conversion\many_subjects_sheet\sub-03\pet	519	50,93	2,887975653	11.15.05	1.017132639884949,1.
D:\BIDS\ONP\BIDS-converter\spreadsheet_conversion\many_subjects_sheet\sub-04\pet	518	19,43	7,555388574	15.45.55	
D:\BIDS\ONP\BIDS-converter\spreadsheet_conversion\many_subjects_sheet\sub-05\pet	469	18,76	7,085	16.12.54	
D:\BIDS\ONP\BIDS-converter\spreadsheet_conversion\many_subjects_sheet\sub-06\pet	496	82,43	1,705282058	12.00.01	
D:\BIDS\ONP\BIDS-converter\spreadsheet_conversion\many_subjects_sheet\sub-07\pet	380	25,57	4,211654282	10.15.05	
D:\BIDS\ONP\BIDS-converter\spreadsheet_conversion\many_subjects_sheet\sub-08\pet	279	42	1,882585714	10.12.36	
D:\BIDS\ONP\BIDS-converter\spreadsheet_conversion\many_subjects_sheet\sub-09\pet	474	29,4	4,569102041	14.25.41	

Blood data

- Convert preformatted xls,xlsx, csv, tsv and bld (PMOD) to blood.tsv and blood.json

A	B	C	D	E	F	G	H
time_manual	manual_whole_blood	manual_plasma	manual_parent_fraction	hplc_recovery_fractions	metadata	time_autosampler	autosampler_whole_blood
5	2,60382574	2,60382574	2,420287093		MetaboliteMethod	0,1	0,180671235
10	1.510.878.819	2,451232644	2,019122324		HPLC	0,2	0,123157239
20	1.250.953.944	2,728974882	1,333615366		MetaboliteRecoveryCorrectionApplied	0,3	0,138270829
40	1.302.358.278	2,930876454	0,713756301		FALSE	0,4	1,558636644
60	1.363.570.232	3,29203593	0,409445452		DispersionCorrected	0,5	12,60661968
					FALSE	0,6	21,64232847
						0,7	28,26041243
						0,8	25,73260612
						0,9	19,43028869
						1	15,76389758
						5	2,60382574
						5,1	1,412707529
						5,2	1,550802421
						5,3	1,983530331
						5,4	1,282172444
						-	-

Blood data

time	plasma_radioactivity	whole_blood_radioactivity	metabolite_parent_fraction	metabolite_polar_fraction
------	----------------------	---------------------------	----------------------------	---------------------------

0	0	0	1	0
---	---	---	---	---

145	43.31	33.79	0.5749	0.1336
-----	-------	-------	--------	--------

292	48.96	37.42	0.3149	0.2746
-----	-------	-------	--------	--------

602	39.84	32.05	0.1469	0.3548
-----	-------	-------	--------	--------

1248	37.38	31.52	0.073	0.444
------	-------	-------	-------	-------

1785	36.40	28.83	0.078	0.429
------	-------	-------	-------	-------

2390	33.13	26.32	0.061	0.453
------	-------	-------	-------	-------

3059	30.83	25.22	0.049	0.473
------	-------	-------	-------	-------

4196	27.28	21.98	0.036	0.503
------	-------	-------	-------	-------

5407	22.70	19.49	0.032	0.523
------	-------	-------	-------	-------

7193	19.71	15.70	0.02	0.559
------	-------	-------	------	-------

```
{
  "PlasmaAvail": true,
  "WholeBloodAvail": true,
  "MetaboliteAvail": true,
  "MetaboliteMethod": "HPLC",
  "MetaboliteRecoveryCorrectionApplied": false,
  "DispersionCorrected": false,
  "time": {
    "Description": "Time in relation to time zero defined by the _pet.json",
    "Units": "s"
  },
  "plasma_radioactivity": {
    "Description": "Radioactivity in plasma samples. Measured using COBRA counter.",
    "Units": "kBq/mL",
    ....
  }
}
```