



***PET-BIDS Tool  
Demonstrations  
and Tutorials***

***Quantification***

**Granville J Matheson**



# Image analysis & Preprocessing



# Quantification

## Considerations

### Study design


- Data and quantification strategy
  - Research question
  - Biological considerations
  - Radiotracer properties

### bloodstream

- Blood
  - Combining blood data
  - Modelling blood data

### kinfitr

- Pharmacokinetic modelling
  - Blood-TAC delay
  - Noise structure (i.e. weights)
  - Model selection
  - Linearised models:  $t^*$  choice



# Blood processing using *bloodstream*

# BIDS Blood Data

	sub-PS26_ses-baselinebrain_rec-DynTOF_pet.json	
	sub-PS26_ses-baselinebrain_rec-DynTOF_pet.nii.gz	
	sub-PS26_ses-baselinebrain_recording-autosampler_blood.json	
	sub-PS26_ses-baselinebrain_recording-autosampler_blood.tsv	
	sub-PS26_ses-baselinebrain_recording-manual_blood.json	
	sub-PS26_ses-baselinebrain_recording-manual_blood.tsv	

PET

Blood


# BIDS Blood Data

Autosampler

Manual

Manual

Manual



# Blood Processing: The Problem

What we want:

- **Frequency:** High temporal granularity
- **Specificity:** Measurements of what we want
- **Accuracy:** Little to no measurement error

In reality

- We *can* measure with high frequency, but not what we really want
- We *can* measure what we want - but with low frequency
- Errors originate from all of the multiple sources

# BIDS Blood Data

## Goals


1. Combine
2. Improve

Autosampler


Manual

Manual

Manual







Manual  
sampling

## Autosampler





### Specification

**Detectors** Well-type NaI(Tl) crystals coupled to a high gain PMT.  
**Analogue** 4096 channels 180° CA/Sec. 11.5MeV mono. automatic dead time correction

Gamma counter

## Whole blood radioactivity





Manual  
sampling




Plasma  
radioactivity

Bound to  
blood cells



**Specification**


**Detectors** Well-type NaI(Tl) crystals coupled to a high gain PMT.  
**Analogue** 4096 channels 500A, 650A, 1150A, 1500A, 1800A, 2048A, 4096A, 8192A, 16384A, 32768A, 65536A, 131072A, 262144A, 524288A, 1048576A, 2097152A, 4194304A, 8388608A, 16777216A, 33554432A, 67108864A, 134217728A, 268435456A, 536870912A, 1073741824A, 2147483648A, 4294967296A, 8589934592A, 17179869184A, 34359738368A, 68719476736A, 137438953472A, 274877906944A, 549755813888A, 1099511627776A, 2199023255552A, 4398046511104A, 8796093022208A, 17592186044416A, 35184372088832A, 70368744177664A, 140737488355328A, 281474976710656A, 562949953421312A, 1125899906842624A, 2251799813685248A, 4503599627370496A, 9007199254740992A, 18014398509481984A, 36028797018963968A, 72057594037927936A, 144115188075855872A, 288230376151711744A, 576460752303423488A, 1152921504606846976A, 2305843009213693952A, 4611686018427387904A, 9223372036854775808A, 18446744073709551616A, 36893488147419103232A, 73786976294838206464A, 147573952589676412928A, 295147905179352825856A, 590295810358705651712A, 1180591620717411303424A, 2361183241434822606848A, 4722366482869645213696A, 9444732965739290427392A, 18889465931478580854784A, 37778931862957161709568A, 75557863725914323419136A, 151115727451828646838272A, 302231454903657293676544A, 604462909807314587353088A, 1208925819614629174706176A, 2417851639229258349412352A, 4835703278458516698824704A, 9671406556917033397649408A, 19342813113834066795298816A, 38685626227668133590597632A, 77371252455336267181195264A, 154742504910672534362390528A, 309485009821345068724781056A, 618970019642690137449562112A, 1237940039285380274899124224A, 2475880078570760549798248448A, 4951760157141521099596496896A, 9903520314283042199192993792A, 19807040628566084398385987584A, 39614081257132168796771975168A, 79228162514264337593543950336A, 158456325028528675187087900672A, 316912650057057350374175801344A, 633825300114114700748351602688A, 1267650600228229401496703205376A, 2535301200456458802993406410752A, 5070602400912917605986812821504A, 10141204801825835211973625643008A, 20282409603651670423947251286016A, 40564819207303340847894502572032A, 81129638414606681695789005144064A, 162259276829213363391578010288128A, 324518553658426726783156020576256A, 649037107316853453566312041152512A, 1298074214633706907132624082305024A, 2596148429267413814265248164610048A, 5192296858534827628530496329220096A, 10384593717069655257060992658440192A, 20769187434139310514121985316880384A, 41538374868278621028243970633760768A, 83076749736557242056487941267521536A, 166153499473114484112975882535043072A, 332306998946228968225951765070086144A, 664613997892457936451903530140172288A, 1329227995784915872903807060280344576A, 2658455991569831745807614120560689152A, 5316911983139663491615228241121378304A, 10633823966279326983230456482242756608A, 21267647932558653966460912964485513216A, 42535295865117307932921825928971026432A, 85070591730234615865843651857942052864A, 170141183460469231731687303715884105728A, 340282366920938463463374607431768211456A, 680564733841876926926749214863536422912A, 1361129467683753853853498429727072845824A, 2722258935367507707706996859454145691648A, 5444517870735015415413993718908291383296A, 10889035741470030830827987437816582766592A, 21778071482940061661655974875633165533184A, 43556142965880123323311949751266331066368A, 87112285931760246646623899502532662132736A, 174224571863520493293247799005065324265472A, 348449143727040986586495598010130648530944A, 696898287454081973172991196020261297061888A, 1393796574908163946345982392040522594123776A, 2787593149816327892691964784081045188247552A, 5575186299632655785383929568162090376495104A, 11150372599265311570767859136324180752990208A, 22300745198530623141535718272648361505980416A, 44601490397061246283071436545296723011960832A, 89202980794122492566142873090593446023921664A, 178405961588244985132285746181186892047843328A, 356811923176489970264571492362373784095686656A, 713623846352979940529142984724747568191373312A, 1427247692705959881058285969449495136382746624A, 2854495385411919762116571938898990272765493248A, 5708990770823839524233143877797980545530986496A, 11417981541647679048466287755595961091061972992A, 22835963083295358096932575511191922182123945984A, 45671926166590716193865151022383844364247891968A, 91343852333181432387730302044767688728495783936A, 182687704666362864775460604089535377456991567872A, 365375409332725729550921208179070754913983135744A, 730750818665451459101842416358141509827966271488A, 1461501637330902918203684832716283019655932542976A, 2923003274661805836407369665432566039311865085952A, 5846006549323611672814739330865132078623730171904A, 11692013098647223345629478661730264157247460343808A, 23384026197294446691258957323460528314494920687616A, 46768052394588893382517914646921056628989841375232A, 93536104789177786765035829293842113257979682750464A, 187072209578355573530071658587684226515959365500928A, 374144419156711147060143317175368453031918731001856A, 748288838313422294120286634350736906063837462003712A, 1496577676626844588240573268701473812127674924007424A, 2993155353253689176481146537402947624255349848014848A, 5986310706507378352962293074805895248510699696029696A, 11972621413014756705924586149611790497021399392059392A, 23945242826029513411849172299223580994042798784118784A, 47890485652059026823698344598447161988085597568237568A, 95780971304118053647396689196894323976171195136475136A, 191561942608236107294793378393788647952342390272950272A, 383123885216472214589586756787577295904684780545900544A, 766247770432944429179173513575154591809369561091801088A, 1532495540865888858358347027150309183618739122183602176A, 3064991081731777716716694054300618367237478244367204352A, 6129982163463555433433388108601236734474956488734408704A, 12259964326927110866866776217202473468949912977468817408A, 24519928653854221733733552434404946937899825954937634816A, 49039857307708443467467104868809893875799651909875269632A, 98079714615416886934934209737619787751599303819750539264A, 196159429230833773869868419475239575503198607639501078528A, 392318858461667547739736838950479151006397215279002157056A, 784637716923335095479473677900958302012794430558004314112A, 1569275433846670190958947355801916604025588861116008628224A, 3138550867693340381917894711603833208051177722232017256448A, 6277101735386680763835789423207666416102355444464034512896A, 12554203470773361527671578846415332832204710888928069025792A, 25108406941546723055343157692830665664409421777856138051584A, 50216813883093446110686315385661331328818843555712276103168A, 100433627766186892221372630771322662657637687111424552206336A, 200867255532373784442745261542645325315275374222849104412672A, 401734511064747568885490523085290650630550748445698208825344A, 803469022129495137770981046170581301261101496891396417650688A, 1606938044258990275541962092341162602522202993782792835301376A, 3213876088517980551083924184682325205044405987565585670602752A, 6427752177035961102167848369364650410088811975131171341205504A, 12855504354071922204335696738729300820177623950262342682411008A, 25711008708143844408671393477458601640355247900524685364822016A, 51422017416287688817342786954917203280710495801049370729644032A, 102844034832575377634685573909834406561420991602098741459288064A, 205688069665150755269371147819668813122841983204197482918576128A, 411376139330301510538742295639337626245683966408394965837152256A, 822752278660603021077484591278675252491367932816789931674304512A, 1645504557321206042154969182557350504982735865633579863348609024A, 3291009114642412084309938365114701009965471731267159726697218048A, 6582018229284824168619876730229402019930943462534319453394436096A, 13164036458569648337239753460458804039861886925068638906788872192A, 26328072917139296674479506920917608079723773850137277813577744384A, 52656145834278593348959013841835216159447547700274555627155488768A, 105312291668557186697918027683670432318895095400549111254310977536A, 210624583337114373395836055367340864637790190801098222508621955072A, 421249166674228746791672110734681729275580381602196445017243910144A, 842498333348457493583344221469363458551160763204392890034487820288A, 1684996666696914987166688442938726917102321526408785780068975640576A, 3369993333393829974333376885877453834204643052817571560137951281152A, 6739986666787659948666753771754907668409286105635143120275902562304A, 13479973333575319897333507543509815336818572211270286240551805124608A, 26959946667150639794667015087019630673637144422540572481103610249216A, 53919893334301279589334030174039261347274288845081144962207220498432A, 107839786668602559178668060348078522694548577690162289924414440996864A, 215679573337205118357336120696157045389097155380324579848828881993728A, 431359146674410236714672241392314090778194310760649159697657763987456A, 862718293348820473429344482784628181556388621521298319395315527974912A, 1725436586697640946858688965569256363112777243042596638790631055949824A, 3450873173395281893717377931138512726225554486085193277581262111899648A, 6901746346790563787434755862277025452451108972170386555162524223799296A, 13803492693581127574869511724554050904902217944340773110325048447598592A, 27606985387162255149739023449108101809804435888681546220650096895197184A, 55213970774324510299478046898216203619608871777363092441300193790394368A, 110427941548649020598956093796432407239217743554726184882600387580788736A, 220855883097298041197912187592864814478435487109452369765200775161577472A, 441711766194596082395824375185729628956870974218904739530401550323154944A, 883423532389192164791648750371459257913741948437809479060803100646309888A, 1766847064778384329583297500742918515827483896875618958121606201292619776A, 3533694129556768659166595001485837031654967793751237916243212402585239552A, 7067388259113537318333190002971674063309935587502475832486424805170479104A, 14134776518227074636666380005943348126619871175004951664972849610340958208A, 28269553036454149273332760011886696253239742350009903329945699220681916416A, 56539106072908298546665520023773392506479484700019806659891398441363832832A, 113078212145816597093331040047546785012958969400039613319782796882727665664A, 226156424291633194186662080095093570025917938800079226639565593765455331328A, 452312848583266388373324160190187140051835877600158453279131187530910662656A, 904625697166532776746648320380374280103671755200316906558262375061821325312A, 1809251394333065553493296640760748560207343510400633813116524750123642650624A, 3618502788666131106986593281521497120414687020801267626233049500247285301248A, 7237005577332262213973186563042994240829374041602535252466099000494570602496A, 14474011154664524427946373126085988481658748083205070504932198000989141204992A, 28948022309329048855892746252171976963317496166410141009864396001978282409984A, 57896044618658097711785492504343953926634992332820282019728792003956564819968A, 115792089237316195423570985008687907853269984665640564039457584007913129639936A, 231584178474632390847141970017375815706539969331281128078915168015826259279872A, 463168356949264781694283940034751631413079938662562256157830336031652518559744A, 926336713898529563388567880069503262826159877325124512315660672063305037119488A, 1852673427797059126777135760139006525652319754650249024631321344126610074238976A, 3705346855594118253554271520278013051304639509300498049262642688253220148477952A, 7410693711188236507108543040556026102609279018600996098525285376506440296955904A, 14821387422376473014217086081112052205218558037201992197050570753012880593911808A, 29642774844752946028434172162224104410437116074403984394101141506025761187823616A, 59285549689505892056868344324448208820874232148807968788202283012051522375647232A, 118571099379011784113736688648896417641748464297615937576404566024103044751294464A, 237142198758023568227473377297792835283496928595231875152809132048206089502588928A, 474284397516047136454946754595585670566993857190463750305618264096412179005177856A, 948568795032094272909893509191171341133987714380927500611236528192824358010355712A, 18971375900641885458197870183





Plasma  
parent

Plasma  
metabolites

Bound to  
blood cells





Manual  
sampling




Plasma  
parent

Plasma  
metabolites

Bound to  
blood cells

HPLC




# BIDS Blood Data

Autosampler


Manual

Manual



Manual






# BIDS Blood Data






Whole Blood Radioactivity



# BIDS Blood Data






# BIDS Blood Data











# BIDS Blood Data




# BIDS Blood Data





# BIDS Blood Data




# BIDS Blood Data





Plasma Parent Fraction




# BIDS Blood Data





Plasma Parent Radioactivity (AIF)



# BIDS Blood Data



Plasma Parent Radioactivity (AIF)



## Basic strategy


Linear interpolation of all curves followed by combination

## More advanced strategy

Model some or all of the functions to reduce measurement error

(Subjective) order of importance:

1. Parent Fraction
2. Blood-to-plasma ratio
3. Arterial Input Function
4. (Whole blood)



# bloodstream

## Usage (docker)

```
docker run --user $(id -u):$(id -g) -v /path/to/bids_data/./data/  
mathesong/bloodstream:latest
```

→ simple linear interpolation of all curves

```
docker run --user $(id -u):$(id -g) -v /path/to/bids_data/./data/  
mathesong/bloodstream:latest config_file.json
```

→ modelling of one or more of the curves



# Usage (R)

```
remotes::install_github("mathesong/bloodstream")
```

→ installation

```
bloodstream("/path/to/bids_data/")
```

→ simple linear interpolation of all curves

```
bloodstream("/path/to/bids_data/", "config_file.json")
```

→ modelling of one or more of the curves

# Create a customised bloodstream config file

## Data Subset

Use these options to apply this config to a subset of the data. Values should be separated by semi-colons. All measurements fulfilling all the conditions will be included. Leave options blank for no subsetting is desired, i.e. leaving sub blank implies that all subjects should be included.

sub

ses

rec

task

run

TracerName

## Modelling Choices

Select the modelling approach for each of the blood curves which should be fitted to the data. The default approach for each is simply to apply linear interpolation to the observed data. As a rule of thumb, modelling the parent fraction and the blood-to-plasma ratio are usually a good idea. Modelling the AIF and the whole blood are mostly best left for specific applications. For debugging, I recommend using simple interpolation and inspecting the plots and QC output.

Parent Fraction

Blood-to-Plasma Ratio

Arterial Input Function

Whole Blood

Download

### Parent Fraction Model Selection

There are many options available for modelling the parent fraction. For most tracers, a good default option is the 'Fit Individually: Choose the best-fitting model' option, which will choose the model which fits best on average, and applies that model to all of the data. Hierarchical models (more to come) are best left for experienced users.

#### Parent fraction model

Interpolation

Fit Individually: Choose the best-fitting model

Fit Individually: Hill

Fit Individually: Exponential

Fit Individually: Power

Fit Individually: Sigmoid

Fit Individually: Inverse Gamma

Fit Individually: Gamma

to (min)

Inf

#### GAM dimension of the basis (k)



6

This value must sometimes be reduced when there are too few data points, or increased for extra wiggleness.

# bloodstream

## Produces

- an html file to inspect all the fits
- Interpolated outputs




# bloodstream

<https://github.com/mathesong/bloodstream>

# Kinetic Modelling using *kinfitr*

# Kinetic modelling



# Kinetic modelling using *kinfitr*

Kinfitr is an R package which contains

- Blood processing methods
  - (bloodstream wraps around these functions)
- Invasive and non-invasive kinetic models
  - ~20 different models
- Weighting schemes
  - For accommodating changing noise the TAC
- Blood-to-TAC delay fitting methods
  - For aligning blood and brain data in time
- $t^*$  selection helpers
  - For using linearised kinetic models
- Miscellaneous
  - Diagnostics: Fit plots, residual plots, ease of model comparison
  - Convenience: Unit conversion, decay correction and un-correction
  - BIDS: PET-BIDS parsing functionality

Full tutorials at  
[granvillematheson.com](http://granvillematheson.com)

## Pharmacokinetic Modelling of PET Data in R using kinfitr. Part 1: Theory.

A brief theoretical background.

*a.k.a. Let me tell you why you this is cool*

*Posted by Granville Matheson Friday, February 21, 2020*

## Pharmacokinetic Modelling of PET Data in R using kinfitr. Part 2: Basics and Iteration

Basic usage and how best to iterate.

*a.k.a. Let's take this baby for a test-drive*

*Posted by Granville Matheson Friday, February 21, 2020*

## Pharmacokinetic Modelling of PET Data in R using kinfitr. Part 3: Finding tstar

How to interpret the tstar finder outputs.

*a.k.a. Admire my busy graphs*

*Posted by Granville Matheson Friday, February 21, 2020*

## Pharmacokinetic Modelling of PET Data in R using kinfitr. Part 4: Blood Processing

Storage and Modelling.

*a.k.a. Model all the things*

*Posted by Granville Matheson Friday, February 21, 2020*

# kinfitr : The BIDS App!

- Needs to be able to automatically parse
  - PET-BIDS data (e.g. blood data)
  - PET-BIDS derivative data (i.e. output from other pipelines)
    - Processed blood data from e.g. *bloodstream*
    - Processed TAC data from e.g. *petprep*
- Needs to allow users to define complex chains of commands and allow checking at intermediate steps
- Want to allow users to experiment with optimal model configurations like graphical user interface (GUI) tools

**Still in active development, but I've released an early version, and working on expanding functionality!**



# kinftr : The BIDS App!

Two parts / interfaces

## 1. Region definition

index	name
2	Left-Cerebral-White-Matter
7	Left-Cerebellum-White-Matter
8	Left-Cerebellum-Cortex
10	Left-Thalamus
11	Left-Caudate
12	Left-Putamen
13	Left-Pallidum
16	Brain-Stem
17	Left-Hippocampus
18	Left-Amygdala
24	CSF
26	Left-Accumbens-area
28	Left-VentralDC
31	Left-choroid-plexus
41	Right-Cerebral-White-Matter
46	Right-Cerebellum-White-Matter
47	Right-Cerebellum-Cortex
49	Right-Thalamus
50	Right-Caudate
51	Right-Putamen
52	Right-Pallidum
53	Right-Hippocampus
54	Right-Amygdala
58	Right-Accumbens-area
60	Right-VentralDC
63	Right-choroid-plexus
130	AirCavity
165	Skull
172	Vermis
174	Pons
257	CSF-ExtraCerebral
258	Head-ExtraCerebral
1001	ctx-lh-bankssts
1002	ctx-lh-caudalanteriorcingulate
1003	ctx-lh-caudalmiddlefrontal
1005	ctx-lh-cuneus
1006	ctx-lh-entorhinal
1007	ctx-lh-fusiform
1008	ctx-lh-inferiorparietal
1009	ctx-lh-inferiortemporal
1010	ctx-lh-isthmuscingulate
1011	ctx-lh-lateraloccipital
1012	ctx-lh-lateralorbitofrontal
1013	ctx-lh-lingual
1014	ctx-lh-medialorbitofrontal
1015	ctx-lh-middletemporal
1016	ctx-lh-parahippocampal
1017	ctx-lh-paracentral
1018	ctx-lh-parsopercularis
1019	ctx-lh-parsorbitalis
1020	ctx-lh-parstriangularis
1021	ctx-lh-pericalcarine
1022	ctx-lh-postcentral
1023	ctx-lh-posteriorcingulate
1024	ctx-lh-precentral
1025	ctx-lh-precuneus

Caudate

Putamen

Cingulate

This is considered common to all  
“analyses” of the data

# kinfitr : The BIDS App!

Two parts / interfaces

1. Region definition
2. Modelling

- Separated into analyses, e.g.
  - Region subsets
  - Session subsets
  - Restricted to first 60 minutes
- Each analysis allows up to 3 models with parameter inheritance
- Each analysis consists of several parts, which can be executed individually, and each of which produces its own report
- A sandbox feature will allow experimenting with model configuration to assess appropriateness

# To the demo...


[https://github.com/mathesong/kinfitr\\_app](https://github.com/mathesong/kinfitr_app)

(usage instructions in the README)

# Delay correction


Input Function Time Shift Optimization

AIF time shift: -2 minutes



Input Function Time Shift and Model Fitting

AIF time shift: -2 minutes



# Weights

